

MY SURECLICK[®] RESOURCE

A GUIDE TO SUPPORT YOUR
REPATHA[®] JOURNEY



YOUR PARTNER IN CARE, EVERY STEP OF THE WAY

This resource is intended to provide you with an overview and basic information about the SureClick autoinjector and does not replace the information in the **Patient Medication Information** leaflet. Please review the instructions for use included with your Repatha[®] (evolucumab injection) prescription prior to performing your injection.

It is important that you do not try to give yourself the injection unless you have received training from your healthcare provider. If you have further questions, please contact your healthcare provider.

 **Repatha**[®]
evolucumab injection

RepathaREADY[®]
PATIENT SUPPORT PROGRAM


by the Amgen Entrust[®] Patient Support Program*

* Amgen Entrust[®] is our unified patient support program platform, built on the legacy of our branded support programs.

GETTING STARTED WITH YOUR SURECLICK® AUTOINJECTOR

This resource is intended to provide you with an overview and basic information about the SureClick autoinjector and does not replace the information in the **Patient Medication Information** leaflet that is included with your Repatha® prescription.

Please review the instructions for use included with your Repatha® prescription prior to performing your injection. If you have further questions, please contact your healthcare provider.



This resource may help you with questions about using the SureClick autoinjector.



See page 20 to learn how to access the SureClick injection demonstration video!

NEED HELP? 

You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.

 **Visit Repatha.ca**

GETTING
STARTED

INJECT

RESOURCES

FAQ

TRAVEL
GUIDANCE

IMPORTANT
SAFETY INFORMATION

GETTING STARTED AT HOME WITH YOUR SURECLICK® AUTOINJECTOR

When you prepare to use the SureClick autoinjector for the very first time, you may wonder if you're going to do it right.

You may have questions or feel uncertain about the process. Other patients feel the same. You are not alone.

There's a lot to learn, and you want to get it right. We're here to support you!

Your Repatha® injection can become a familiar part of your regular routine – one injection, every two weeks.



With this resource and the online tools, **you can do this.**
Now, let's get started!

Carefully review **My SureClick Resource**, which offers additional helpful information just for you and your caregiver, including:

1. Key injection steps
2. Helpful tips
3. Resources and tools
4. Answers to frequently asked questions

Please see important safety information on page 25. Refer also to the full Patient Medication Information leaflet, which is provided with your medication.

GETTING
STARTED

INJECT

RESOURCES

FAQ

TRAVEL
GUIDANCE


IMPORTANT
SAFETY INFORMATION


FIND A QUIET, CALM PLACE TO INJECT

Start off right and don't rush

Find a quiet and calm place in your home, where pets or other people can't distract or interrupt you.

Consider a room where you can close the door. Sit in a comfortable chair so you feel more relaxed.



NEED HELP?  You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.

 [Visit Repatha.ca](https://www.repatha.ca)

INJECT

In this section, we will discuss the steps for injecting the medicine.

Please read the step-by-step instructions in the **Patient Medication Information** leaflet prior to starting the injection.



Please see important safety information on page 25. Refer also to the full Patient Medication Information leaflet, which is provided with your medication.

GETTING
STARTED

INJECT

RESOURCES

FAQ

TRAVEL
GUIDANCE

IMPORTANT
SAFETY INFORMATION


GET TO KNOW YOUR SURECLICK AUTOINJECTOR

An introduction: Important things to know about your SureClick® autoinjector – Before Use

- SureClick autoinjector is a single-use autoinjector prefilled with your medicine.

⚠ Important: Before you use the SureClick autoinjector, refer to the **Patient Medication Information** leaflet that comes with your Repatha medication for **important information** about **using and storing** the SureClick autoinjector, along with **step-by-step instructions**.

Do not try to inject Repatha until you have been shown the right way by a healthcare professional, such as a doctor, nurse or pharmacist.




⚠ Do not use the autoinjector if the medicine is cloudy or discoloured or contains large lumps, flakes or particles.

⚠ Do not use the autoinjector if the expiration date printed on the label has passed

⚠ Do not remove orange cap until ready to inject. Inject within **5 minutes** of removing orange cap.



- If the orange cap is off for more than 5 minutes, this can dry out the medicine and may clog the needle.

NEED HELP?  You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.


 [Visit Repatha.ca](http://Repatha.ca)

GET TO KNOW YOUR SURECLICK AUTOINJECTOR

An introduction: Important things to know about your SureClick® autoinjector – After Use



TURNING YELLOW





COMPLETE ✓

WINDOW will turn **yellow** during the injection (injection complete)

Make sure the medicine in the window is clear and colourless to slightly yellow


The **YELLOW SAFETY GUARD** contains the **needle**. The guard must be pushed down on the skin until it stops moving in order to release the needle when the **GREY START BUTTON** is pressed. **Do not** touch the grey button until you are ready to inject.

YELLOW SAFETY GUARD IS **NOT** PUSHED ALL THE WAY DOWN



YELLOW SAFETY GUARD **IS** PUSHED ALL THE WAY DOWN


 **Do not** put fingers into the yellow safety guard.

NEED HELP?  You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.

 Visit Repatha.ca

PREPARE | READY THE AUTOINJECTOR


Prepare the autoinjector by bringing it to room temperature naturally



Remove the autoinjector from the carton

- Carefully lift the autoinjector straight up out of the carton.
- Before you inject, always check the label of the autoinjector to make sure you have the correct medicine and the correct dose of Repatha.
- Put the original carton with any unused autoinjectors back in the refrigerator.
- Gather all materials needed for your injection: alcohol wipes, cotton balls or gauze pads, adhesive bandages and a sharps disposal container.

 Keep the autoinjector and all medicines out of the reach and sight of children.




Wait at least 30 minutes for the autoinjector to reach room temperature before injecting.

This is important for administering the entire dose and helps minimize discomfort. Repatha may take longer to inject if it has not reached room temperature. Let it come to room temperature naturally.




Pull the orange cap straight off only when you are ready to inject.

-  **Do not** leave the orange cap off for more than 5 minutes
- This can dry out the medicine and may clog the needle when you take the injection.


PREPARE | READY THE AUTOINJECTOR




 **Do not** shake the autoinjector.


 Do not store the autoinjector in extreme heat or cold. Store in a refrigerator at 2°C to 8°C in the original carton.

- When removed from the refrigerator, Repatha should be kept at controlled room temperature (up to 25°C) in the original carton and must be used **within 30 days**.


- Throw away REPATHA that has been stored at room temperature for more than 30 days.


 **Do not** leave the autoinjector in the glove box or trunk of your car.




 Protect Repatha from direct light and do not expose to temperatures above 25°C.



 **Do not** try to warm the autoinjector using a heat source such as hot water or microwave.



 **Do not** freeze or use the autoinjector if it has been frozen.

PREPARE | SELECT AND CLEAN THE INJECTION SITE


Select the injection site for you or your caregiver

- Begin by gathering materials in a quiet, calm, and well-lit area. Choose a location to ensure you aren't interrupted. Place supplies on a flat surface.

- **Only use these injection sites**

- Thigh
- Stomach (except for a 5-centimetre [2-inch] area around your belly button)
- Outer area of the upper arm (only if someone else is giving you the injection).

- Wash hands before selecting and cleaning the injection site.



- Clean your injection site with an alcohol wipe and let your skin dry

⚠ Do not touch this area again before injecting.

- Choose a different injection site each time you give yourself an injection. If you want to use the same injection site, make sure it is not the same spot you used for the last injection.
- **⚠ Do not** inject into areas where skin is tender, bruised, red, or hard. Avoid injecting into areas with scars or stretch marks.
- Avoid injecting directly into raised, thick, red, or scaly skin patch or lesion


NEED HELP?  You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.

 [Visit Repatha.ca](https://www.Repatha.ca)

GET READY | REMOVE THE CAP FROM THE AUTOINJECTOR

Pull the orange cap straight off only when you are ready to inject

If you have questions, seek help by referencing the **Patient Medication Information** leaflet or calling the RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST. **Or visit Repatha.ca**




- **Do not** remove the orange cap from the autoinjector until you are ready to inject. If you are unable to inject, please contact your healthcare provider.

- ⚠ **Do not** leave the orange cap off for more than 5 minutes. This can dry out the medicine and may clog the needle when you give the injection.

- ⚠ **Do not** twist, bend, or wiggle the orange cap.
- ⚠ **Do not** put the orange cap back onto the autoinjector once it has been removed.
- ⚠ **Do not** put fingers into the yellow safety guard.


GET READY | PINCH OR STRETCH THE SKIN TO CREATE AN AREA ABOUT 5 CENTIMETRES (2 INCHES) WIDE

STRETCH METHOD 2 inches (5 cm)




Stretch the skin firmly by moving your thumb and fingers in opposite directions, creating an area about **5 centimetres (2 inches) wide**.

PINCH METHOD 2 inches (5 cm)




Pinch the skin firmly between your thumb and fingers, creating an area about **5 centimetres (2 inches) wide**.

 It is important to keep skin stretched or pinched throughout the injection.


NEED HELP?  You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.

 [Visit Repatha.ca](https://www.Repatha.ca)


INJECT | PLACE THE SAFETY GUARD ON THE SKIN AT A 90-DEGREE ANGLE



Keep holding the stretched or pinched skin. With the orange cap off, **put** the yellow safety guard on your skin at 90 degrees.



Be sure you can see the entire medicine window and that your fingers don't cover the window.




⚠ Do not apply the autoinjector to the skin at an angle. You may not be able to make full contact with the skin, and that can prevent the autoinjector from working properly.

- ⚠ Do not** touch the grey start button yet.
- ⚠ Do not** change positions after the injection has begun.

INJECT | PUSH AUTOINJECTOR FIRMLY ONTO SKIN AND HOLD


To start the injection, **PUSH** the autoinjector firmly onto the skin until it stops moving and the yellow safety guard is no longer visible



PUSH straight down and **HOLD**

- Firmly **push** autoinjector down onto the skin until it stops moving.
- Make sure the yellow safety guard is fully pushed down before the next step.
- **Hold** until injection completes.


- ⚠ The **GREY START BUTTON** will work only after the yellow safety guard has been pushed all the way down.
- ⚠ Keep pushing straight down throughout the entire injection.

NEED HELP?  You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.

 [Visit Repatha.ca](https://www.repatha.ca)

INJECT | PRESS GREY START BUTTON TO BEGIN THE INJECTION

Next press the grey start button when you are ready to start the injection




PRESS grey start button until you hear a click.

- **Keep holding** the autoinjector down with the same pressure throughout the injection.
- **GREY START BUTTON** will work only after the yellow safety guard has been pushed all the way down.

 These steps need to be completed in this sequence for a successful injection.

INJECT | WATCH AND CONFIRM MEDICINE WINDOW TURNS **YELLOW**


Finally, watch and confirm the medicine window turns **yellow** to indicate the successful completion of the injection



WATCH AND CONFIRM

window is fully **yellow**


- Keep **pushing** the autoinjector down on the skin. Then **lift** your thumb while still holding the autoinjector on your skin. The injection may take up to 15 seconds to complete.
 - Window turns from clear to yellow when the injection is done. You may hear a second click.
- Window turns from clear to **yellow** when the injection is done. You may hear a **second click**. (**Note:** after you remove the autoinjector from your skin, the needle will be automatically covered.)

NEED HELP?  You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.

 [Visit Repatha.ca](https://www.repatha.ca)

STORAGE AND DISPOSAL

Refrigeration guidelines:



⚠ **Do not** freeze or use the autoinjector if it has been frozen.

Do not store in extreme heat or cold.

- When you receive the autoinjector, be sure to place it directly into the refrigerator.
- Keep the autoinjector in the refrigerator at 2° to 8°C.
- Keep the autoinjector in original carton to protect from light.


Disposal guidelines:

- Put the used autoinjector and orange cap in a sharps disposal container right away after use.

⚠ **Do not** throw away (dispose of) the autoinjector or orange cap in your household trash.

⚠ **Do not** recap the autoinjector.

⚠ **Do not** put fingers into the yellow safety guard.



CONGRATULATIONS

on reviewing the injection process!

Continue to practice the steps **PUSH, PRESS, WATCH AND CONFIRM** to help build confidence. Start to make these steps a part of your Repatha® treatment routine, and be sure to contact your care team if you have any questions or concerns.



NEED HELP? 

You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.



Visit [Repatha.ca](https://www.Repatha.ca)

RESOURCES

In this section, we will discuss a variety of helpful resources available to you that are just a click or phone call away.



Please see important safety information on page 25. Refer also to the full Patient Medication Information leaflet, which is provided with your medication.

GETTING
STARTED

INJECT

RESOURCES

FAQ

TRAVEL
GUIDANCE

IMPORTANT
SAFETY INFORMATION

SURECLICK AUTOINJECTOR PATIENT SUPPORT IS AVAILABLE

Live support



RepathaREADY® by the Amgen Entrust® Patient Support Program is ready to help.

Call 1-888-Repatha (1-888-737-2842), Monday to Friday, 8:00 am to 8:00 pm EST.

Repatha® SureClick® autoinjector demonstration video



Hold your smartphone camera over the QR code.


The Repatha SureClick autoinjector demonstration video can also be viewed at **Repatha.ca**

Additional resources



Visit **Repatha.ca**

To enter, you will need the SureClick autoinjector DIN: 02446057.

NEED HELP?  You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.

 **Visit Repatha.ca**

TEMPERATURE CONTROL AND STORAGE:

What is the best way to store my SureClick® autoinjector? Store in the refrigerator (2°C to 8°C) in its original carton.

⚠ Do not freeze.

I left my autoinjector outside of the refrigerator. Is this okay? After removing from the refrigerator, the autoinjector should be kept at controlled room temperature (up to 25°C) in the original carton and must be used within 30 days.

- **Do not** return the autoinjector to refrigerator once it has reached room temperature.

Room temperature:
Up to 25°C Up to 30 days




ROOM TEMPERATURE GUIDELINES:

Do I really have to wait 30 minutes for the autoinjector to reach room temperature before injecting?

Yes. This is important for administering the entire dose and helps minimize discomfort. Repatha® may take longer to inject if it has not reached room temperature. **Do not** heat the autoinjector. Let it come to room temperature naturally.

- **Do not** try to warm the autoinjector by using a heat source such as hot water or microwave.

Why do I have to use the autoinjector within 5 minutes after I remove the cap? If you wait longer than 5 minutes, this can dry out the medicine and may clog the needle when you give the injection.



Please see important safety information on page 25. Refer also to the full Patient Medication Information leaflet, which is provided with your medication.

GETTING STARTED

INJECT

RESOURCES

FAQ

TRAVEL GUIDANCE

IMPORTANT SAFETY INFORMATION

FAQ

INJECTION SITE QUESTIONS:

Why do I need to stretch or pinch the skin? You need to **stretch** or **pinch** the skin to create an area 5 centimetres (2 inches) wide **to maintain a firm base** throughout the entire injection.

Can I move the autoinjector around on my skin? It is okay to move the autoinjector around on the injection site as long as you **do not** press the grey start button. However, if you press the grey start button and the yellow safety guard is pushed into the autoinjector, the injection will begin.


⚠ Do not relocate the autoinjector **after** you **push** the grey **start** button, as you may not receive the full dose of medicine.

START BUTTON QUESTIONS:

Do I press and release the grey start button or do I press and hold down the grey start button during the injection? When you are ready to inject, **press** the grey start button. You will hear a click. Keep pushing the autoinjector down on the skin. Then **lift** your thumb while still holding the autoinjector on your skin. The injection may take up to 15 seconds to complete.

What will happen if I press the grey start button before I am ready to do the injection on my skin? The injection will only happen when the grey start button is pressed and the autoinjector is pushed into the injection site. Release the grey start button and push the autoinjector down until the yellow safety guard stops moving. Then, you can push the grey start button again.

What if I am unable to push down the grey start button? Check to see that the yellow safety guard is no longer visible. If the yellow safety guard is not depressed fully, the grey start button may not work. You will hear a **CLICK** when the grey start button has been activated correctly.



TRAVEL GUIDANCE

GENERAL TRAVEL GUIDANCE

Yes, you can travel. You don't have to miss any of your treatments. In fact, it's important that you don't miss your injection when travelling.

- Always keep your medicine with you to ensure it is kept within the safe temperature zone and is not broken.

⚠ Do not store the autoinjector in extreme heat or cold. For example, avoid storing the autoinjector in your vehicle's glove box or trunk.

⚠ Do not shake the autoinjector.

⚠ Do not use the autoinjector if it has been dropped on a hard surface. Part of the autoinjector may be broken even if you cannot see the break. Use a new autoinjector, and call **1-888-Repatha (1-888-737-2842)**.

SHORT TRIPS

- If removed from the refrigerator, the autoinjector should be kept at controlled room temperature (up to 25°C) in the original carton and must be used within 30 days.

⚠ Remember that once the autoinjector reaches room temperature, you should not put it back in the refrigerator.

Important:
Keep the autoinjector in the original carton.



LONGER TRIPS

- Reach out to your pharmacist and/or doctor to discuss other options such as filling your prescription at your desired destination.



Please see important safety information on page 25. Refer also to the full Patient Medication Information leaflet, which is provided with your medication.

GETTING
STARTED

INJECT

RESOURCES

FAQ

TRAVEL
GUIDANCE

IMPORTANT
SAFETY INFORMATION

TRAVEL GUIDANCE


CAR AND AIR TRAVEL GUIDANCE

Once the autoinjector reaches room temperature, do not put it back in the refrigerator.

- Throw away REPATHA that has been stored at room temperature for more than 30 days.

Travelling by car:

 **Do not** leave the autoinjector in your car if exposed to extreme heat, extreme cold, or direct sunlight.




 **Do not** store the autoinjector in your vehicle's glove box or trunk.

Travelling by airplane:

Before you travel:

- Contact your airline about the carry-on policy for needles, syringes, or injection devices. Security regulations may change.
- Research the injection device requirements of other countries you may visit.
- Bring a copy of your prescription or a doctor's note to show airport security.



Airport security considerations:

- Notify security screeners that you are carrying a sharps container.
- Consider requesting a "hand check" of your medicine kit to avoid any delays.
- Be prepared to display, handle, and repack your own medicine and supplies in order to prevent contamination or damage.



For information on airline travel with Repatha®, please refer to: Canadian Air Transport Security Authority (CATSA) www.catsa-acsta.gc.ca/

- If you can carry-on Repatha, ensure the **prescription (Rx) label** is attached.
- If you need extra help or have concerns about airport security, please call **1-888-294-2202**.
- For flights originating **outside of Canada**, check with your air carrier or travel agent for any restrictions.

NEED HELP? 

You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.



Visit Repatha.ca

IMPORTANT INFORMATION ABOUT REPATHA®

Do not use Repatha if you have ever had an allergic reaction to Repatha or any of the ingredients in Repatha.

To help avoid side effects and ensure proper use, talk to your healthcare professional before you take Repatha. Talk about any health conditions or problems you may have.

Before starting Repatha, you or your child should be on a diet to lower cholesterol. If your doctor has prescribed Repatha along with a statin or other cholesterol lowering medicine for you or your child, follow your doctor's instructions on how to take these medicines together. In this case, please read the dosage instructions in the package leaflet of the other medicines.

Children and adolescents

The use of Repatha has been studied in children 10 years of age and older being treated for heterozygous or homozygous familial hypercholesterolemia.

The use of Repatha has not been studied in children under 10 years of age or in children with other types of hypercholesterolemia.

Other medicines and Repatha

Tell your doctor or pharmacist about all the medications you are taking or have recently taken, including medicines obtained without a prescription such as drugs, vitamins, minerals, natural supplements or alternative medicines.

Pregnancy and breastfeeding

Repatha has not been tested in pregnant women. It is not known if Repatha will harm your unborn baby. If you are trying to get pregnant or become pregnant while taking Repatha:

- Inform your doctor.
- If you are also taking a statin along with Repatha, stop taking Repatha and read the package leaflet of the statin that you are taking with Repatha.

It is not known whether Repatha is found in breast milk. It is important to tell your doctor if you are breastfeeding or plan to do so. Your doctor will then help you decide whether to stop breastfeeding, or whether to stop taking Repatha, considering the benefit of breastfeeding to the baby and the benefit of Repatha to the mother.

Please see important safety information on page 25. Refer also to the full Patient Medication Information leaflet, which is provided with your medication.

GETTING
STARTED

INJECT

RESOURCES

FAQ

TRAVEL
GUIDANCE

IMPORTANT
SAFETY INFORMATION

WHAT IS REPATHA® AND WHAT IS IT USED FOR?

Repatha is a medicine used to lower levels of cholesterol.

Why have I been prescribed Repatha?

Repatha is used:

- along with diet and in combination with usual therapy, including cholesterol-lowering medications called statins, to reduce the risk of heart attack, stroke and certain heart procedures (to restore blood flow to the heart) in adult patients who have cardiovascular disease by further lowering low-density lipoprotein (LDL) cholesterol.
- along with diet, alone or together with other cholesterol-lowering

therapies, in adults with primary hyperlipidemia (high LDL cholesterol level in your blood) to reduce LDL cholesterol.

- along with diet and other cholesterol-lowering therapies, to reduce LDL cholesterol in people 10 years and older who need additional lowering of the high LDL cholesterol in their blood because of a condition that runs in their family (heterozygous familial hypercholesterolemia [HeFH]).
- along with diet and other cholesterol-lowering therapies, to reduce LDL cholesterol in people 10 years and older who need

additional lowering of the high LDL cholesterol in their blood because of a condition that runs in their family (homozygous familial hypercholesterolemia [HoFH]).

NEED HELP? 

You can call RepathaREADY® by the Amgen Entrust® Patient Support Program for assistance at **1-888-Repatha (1-888-737-2842)**, Monday to Friday, 8:00 am to 8:00 pm EST.



Visit [Repatha.ca](https://www.repatha.ca)



RepathaREADY® Patient Support Program is ready to help.

by the Amgen Entrust® Patient Support Program*

**Call 1-888-Repatha (1-888-737-2842),
Monday to Friday, 8:00 am to 8:00 pm EST.**

Visit Repatha.ca

(SureClick autoinjector
DIN: 02446057) or scan
the QR code.



* Amgen Entrust® is our unified patient support program platform, built on the legacy of our branded support programs. © 2025 Amgen Canada Inc. All rights reserved. Repatha®, RepathaREADY® and SureClick® are registered trademarks of Amgen Inc., used with permission. AMGEN ENTRUST® is a registered trademark of Amgen Inc.

